

3M Health Care and Materials Management Microsystems Announce Launch of 3M Attest Monitoring Interface for SPM Software

- Interface Delivers Enhanced Connectivity and Efficiency for Sterile Processing -

3M Health Care and Materials Management Microsystems today announced the launch of 3M™ Attest™ Monitoring Interface for SPM® Software, designed to provide connectivity between Microsystems' SPM® Sterile Processing Information System and the 3M Attest 290 and 290G Auto-readers. This connectivity provides hospitals more streamlined data collection, documentation and reporting capabilities as results from the 3M Attest Rapid Readout biological indicator (BI) tests download directly into Microsystems' software application.

Through this offering, Microsystems provides a new module to its existing software application that includes enhanced sterilization monitoring documentation capabilities and instant access to real-time results. Sterile processing clients will enjoy greater accuracy and simplicity in the documentation process as more information moves electronically to reside in a single, convenient location.

"We are thrilled to have an opportunity to improve the efficiency of the Sterile Processing Department (SPD) by offering instant access to 3M Attest Rapid Readout biological indicator results and real-time documentation of those results," said Tammy Torbert, 3M Sterilization Assurance marketing manager, 3M Health Care. "As scrutiny toward hospital-acquired infections (HAIs) continues to increase, the electronic documentation of critical sterilization monitoring results will allow the SPD to demonstrate sound quality control processes."

Donna McLaughlin, Sterile Processing Manager at Fort Sanders Regional Medical Center, a facility in Knoxville, TN that recently served as a test site for the new Microsystems/3M offering, affirms the importance of quality control processes in her line of work: "We are in the business of providing excellent patient care and that level of care is synonymous with patient safety. My department supports patient safety initiatives by sterilizing surgical instruments, and it is imperative that we have the highest level of monitoring controls in place to do that. Having the 3M Attest Monitoring Interface for SPM Software helps us manage every load monitoring more efficiently; we run a BI in every steam load, and the interface gives us real-time results with less potential for human error."

Microsystems echoes this praise and support of the new connectivity: "The 3M Attest Monitoring Interface for SPM Software represents a significant advancement in our efforts to combine essential information from multiple sources," said Clay Cannady, Microsystems' head of Sales, Marketing and Business Development. "The integration of our companies' proven, data-driven solutions arms sterilization professionals with the critical information they need when assessing sterility for surgical procedures."

The 3M Attest Monitoring Interface for SPM Software is now available. 3M and Microsystems are collaborating on the sales and marketing of this new offering. As such, customers can learn more about this product from their 3M or Microsystems sales representative. Facilities may also contact Microsystems at 1.262.240.9900, Option 4, or sales@mmmicrosystems.com.

About 3M Sterilization Assurance

With over 40 years of providing sterilization solutions, 3M Sterilization Assurance is a trusted choice for

sterilization professionals. 3M Sterilization Assurance, a part of 3M Health Care's infection prevention platform, offers well known products including 3M™ Attest™ Rapid Readout Biological Indicators, 3M™ Comply™ SteriGage™ Chemical Integrators, and 3M™ Steri-Vac™ Sterilization Equipment to the healthcare market.

About 3M Health Care

3M Health Care, one of 3M's six major business segments, provides world-class innovative products and services to help healthcare professionals improve the practice, delivery and outcome of patient care in medical, oral care, drug delivery and health information markets. Learn more at www.3m.com/healthcare.

3M, Attest, Comply, SteriGage and Steri-Vac are trademarks of 3M.

About Materials Management Microsystems

Microsystems provides bar code driven sterile processing information systems designed specifically for healthcare facilities. Microsystems is a leader in this market by delivering a unique combination of sophisticated products, refined services, attentive and personal customer support, and a staff of professionals committed to the success of each client. Learn more at www.mmmicrosystems.com.

SPM is a trademark of Materials Management Microsystems.

Materials Management Microsystems
Janet Becker, 262-240-9900 x900
jbecker@mmmicrosystems.com
or Mary Clemens, 651-733-8806
mcclemens@mmm.com

<https://news.3m.com/2008-12-08-3M-Health-Care-and-Materials-Management-Microsystems-Announce-Launch-of-3M-Attest-Monitoring-Interface-for-SPM-Software>